

Del Medical— Delivering Quality Products

to the imaging profession for more than 75 years.

Quality and cost-effective performance

Consistently ranked high for exceptional performance, Del Medical offers uncompromising radiographic systems that are an excellent economical alternative for hospitals, clinics and private practices looking for ease-of-use, image quality and high patient throughput.

Future-safe and digital-ready

As the field of radiology changes and advances, count on Del Medical to respond with technological innovations that support those advances, and with systems that operate with a low life-cycle cost and offer the highest productivity. Del Medical systems are "DR-ready," giving our customers the option of delivering improved image quality and workflow for a full range of general radiographic imaging applications.

Service beyond the sale

Service is never an afterthought with Del Medical, but an integral part of the purchasing agreement. Your investment in Del Medical equipment is always backed by our company's commitment to maintain the equipment long after installation. Our worldwide dealer network, expertly trained on our product lines, provides prompt service and routine preventive maintenance.

Trust the leaders in radiology solutions

For decades, conventional radiographic systems have been the mainstay of diagnostic imaging at virtually every hospital or medical center. Del Medical has been a leader in this market for years — evolving with the changing needs of the industry, providing the most technologically advanced equipment available. Del Medical is the leader it is today for one simple reason — our customers recognize that we deliver the best return on their investment, system after system.




Accessories and Options



Floor-mounted tubestand



Hand grips



Lateral cassette holder



Made in the USA

ADVANCED RADIOGRAPHIC SYSTEM



A Durable and Easy-To-Operate Platform
For Your Radiographic Future

DEL MEDICAL

Rev. 8/05 #9000-EV650/OTC3

USA: 11550 West King Street • Franklin Park, IL 60131
EUROPE: Via Dell Azalee, 3 • 20090 Buccinasco • Milan, Italy
Phone 800.800.6006 • Fax 847.288.7011

www.delmedical.com

DEL MEDICAL

What do you need from your radiographic system?

Whether it's conventional chest x-ray, trauma imaging or a variety of other general radiographic studies, Del Medical's systems are up to the task. With Del, you are assured of a radiographic room that perfectly fits your hospital's or clinic's workflow and budget. You will also have the freedom to 'scale-up' to new functionality – Del radiographic systems come fully upgradable to CR or DR digital imaging.

The ideal system
for fast, efficient and
operator-friendly workflow.



Rugged design for high volume environments

Built to excel in the most demanding environments and to hospital-quality duty standards, Del Medical's EV 650 is a rugged four-way float top elevating table designed to produce high quality images with effortless patient transfer and positioning.

Effortless head-to-toe imaging

Del's EV 650 four-way float top elevating table provides for smooth, quick and accurate patient positioning. The EV 650 takes radiographic systems to new levels of durability and safety with its 650 lb (295kg) patient weight capacity, special high strength contoured table façade, integral collision sensors and recessed, slanted foot treads that reduce accidental activation. The EV 650 also incorporates the largest, quietest and most powerful elevating assembly in the industry to insure reliable operation, and will provide for years of trouble-free operation.



Designed to accommodate any study

The OTC₃ overhead tube stand extends the versatility of Del Medical's radiographic systems. By providing greater freedom of movement within the room, it facilitates examinations performed on any radiographic study. Further, the OTC₃ has a center mounted single-handle lock release that enables the positioning of the tube assembly for quick positioning in emergency situations.

Del Medical's OTC₃ is designed for minimum effort and ease of positioning, as it moves smoothly along an overhead structure that rides on special bearing assemblies. To assist in procedural set up, digital readouts are provided along with easy-touch switches. Its unique design reduces installation time and costs. The OTC₃ is available in several configurations designed to accommodate your application requirements.

