## **KODAK DryView 8900 Laser Imager**

For SOMATOM VolumeZoom, Volume Access, Sensation 4, Sensation 10, Sensation 16, Sensation Cardiac, Esprit, Balance, Emotion, Emotion Duo, Emotion 6.

The DryView 8900 has remarkable speed and superb image quality, making it an excellent seller worldwide. It renders sharp alphanumeric and continuous-tone images with exceptional spatial resolution, contrast and greyscale.

## Includes

- DryView 8900 Laser Imager prints film sizes (14"x17", 14"x14", 11"x14", 10"x12", 8"x10")
- Six (6) connections: DICOM or Non – DICOM
- Top mounted film sorter six bins (sorts by modality)
- Three film drawers (each drawer holds any size film, 125 sheets/cartridge)

## Features

- Up to 200 sheets of film per hour every image at 650 dpi
- 39-micron laser spot, 14 bit
- Continuous-tone images at 16,384 gray levels
- Connects to all DICOM network modalities and serves 12 simultaneously
  - Kodak's Automatic Image Quality Control maintains film to film consistency
- Increased capability in multi-format imaging without compression or loss of data
  - Large, intuitive, color touch screen

Includes installation, applications and 1 year warranty through KODAK.



Order No.

K114 7107

